ACCU-CHEK, ACCU-CHEK FLEXLINK, ACCU-CHEK LINKASSIST, ACCU-CHEK SPIRIT, ACCU-CHEK SPIRIT COMBO, and COMBO are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

© 2023 Roche Diabetes Care

Made in Mexico

Roche Diabetes Care GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany
www.accu-chek.com

Last update: 2023-05

04537246001/Q


ACCU-CHEK® FlexLink

Roche

Instructions for use

Infusion set

EN DE FR IT ES PT / PTB NL
DA NO SV FI PL CS

ACCU-CHEK®

PDF grossed3 - Jul 13, 2023

Approved according to
SHR SOP 04.07.08 Att. A Vers.5.3


Kat. 03


The Signature will not be printed!

GD/
LAY grossed3 - Jul 13, 2023


PM/
LAB hoffartm - Jul 19, 2023

LAN caspersc - Jul 13, 2023

1

Suitable infusion sites
Geeignete Infusionsstellen
Sites de perfusion appropriés
Siti d'infusione idonei
Zonas de infusión adecuadas
Locais de infusão adequados
Geschikte infusieplaatsen
Passende infusionssteder
Egnede infusionssteder
Lämpliga infusionsställen
Sopivat infuusiokohdat
Odpowiednie miejsca infuzji
Vhodná místa infuze

Instructions for use	4	en
Gebrauchsanweisung	16	de
Instructions d'utilisation	30	fr
Istruzioni per l'uso	44	it
Instrucciones de uso	58	es
Instruções de utilização/uso	72	pt/ptb
Gebruiksaanwijzing	86	nl
Brugsanvisning	100	da
Brugsanvisning	112	no
Bruksanvisning	124	sv
Käyttöohje	136	fi
Instrukcja obsługi	148	pl
Návod k použití	162	cs
	3	

PRECAUTION

Risk of serious harm

The disregard of the instructions for use can prevent the product from working as intended. This can lead to various harms, including hypoglycaemia, hyperglycaemia, injury, or infection.

Carefully read and follow the instructions for use. Also, read the User's Manual of your insulin pump.

Intended Use

The infusion set is intended for the subcutaneous infusion of insulin delivered by an insulin pump.

Indication

For use in the therapy of diabetes mellitus.

Intended User

People with diabetes and healthcare professionals or caregivers who assist people with diabetes in using the product.

Target Group (applicable for)

Infants, children, adolescents, and adults with diabetes.

Contraindications

There is no clinical experience for the use with neonates.

Clinical Benefits

The infusion set enables pump therapy when it is connected to a compatible insulin pump. The infusion set enables connected insulin pumps to fulfil their purpose, which is to provide an intensive insulin therapy.

Limitations

The device is intended for single use only.
The device is intended to be used only with insulin or insulin analogues.

Additional Materials Required

The infusion set is intended to be used together with the following products only:

Insulin pumps

- the Accu-Chek Spirit insulin pump
- the Accu-Chek Spirit Combo insulin pump

Cartridges

- the Accu-Chek Spirit 3.15 ml cartridge system
- the Insuman® Infusat 100 IU/ml solution for injection in a cartridge (3.15 ml, Sanofi Aventis)

Optional Accessory

- the Accu-Chek LinkAssist insertion device

Device Description

Head set

- A** Handling aid
- B** Side wings
- C** Adhesive pad
- D** Protective film
- E** Cannula
- F** Introducer needle
- G** Needle cover
- H** Cannula casing
- I** Disconnect cover

Transfer set

- J** Luer-lock connector
- K** Tubing

en

- L** Connector
- M** Side wings

The infusion set includes a head set and a transfer set. The head set contains a soft cannula and the introducer needle. The pointed and sharp introducer needle inserts the cannula at a 90° angle into the subcutaneous tissue.

The transfer set must be connected to the cannula and through a standard luer-lock connector to the insulin pump.

The head set and the transfer set are sterile, single use products that have been sterilised with ethylene oxide.

Contents of the Pack

There are 2 packages with different contents.

Infusion set: Contains head sets, disconnect covers, and transfer sets, **[REF]** 04626427,

04626435, 04626478, 04626486, 04626494, 04626516, 04626524, 04626559, 06593844, 06593909, 06593917

Infusion set cannula: Contains head sets and disconnect covers, **[REF]** 04626567, 04626575, 06593941

General Safety Information

A  **WARNING** indicates a foreseeable serious hazard.

A  **PRECAUTION** describes a measure you should take to use the product safely and effectively or to prevent damage to the product.

WARNING

Risk of suffocation

This product contains small parts that can be swallowed.

Keep the small parts away from small children

and people who might swallow small parts.

Risk of strangulation

If the tubing of the transfer set is wrapped around the neck, the tubing can strangle you. Pay attention to the tubing. The transfer set is available in different lengths. Use the shortest possible length.

Risk of hyperglycaemia

If the connections in the infusion set are not tight, the insulin delivery can be insufficient. If you feel unwell, if you suspect insufficient insulin delivery, or if you unexpectedly smell insulin, check whether your blood glucose value increases. Also, check the infusion set for leaks.

Risk of hyperglycaemia or infection

Expired products can lead to an insufficient insulin delivery or can be unsterile. Use only products that are within the use by

date. The use by date is printed next to the  symbol on the packaging.

Risk of infection

If the sterile packaging is damaged or already open, the product can be unsterile and lead to infection.

Visually inspect the sterile packaging. If the packaging is damaged or already open, do not use the product. The sterile packaging is marked with the  symbol.

PRECAUTION

Risk of hyperglycaemia, increased pain, or injury

If you insert the cannula into a body part that is not suitable as infusion site, the insulin delivery can be insufficient or increased pain or injury can occur.

Select a suitable infusion site. Ask your

en

7

healthcare professional where you can insert the cannula.

Risk of hyperglycaemia

- Several reasons can lead to an insufficient insulin delivery.

After attaching or replacing the cannula or the transfer set, check your blood glucose level at least once within 1 to 3 hours. If your blood glucose level is unexpectedly high, follow the instructions in these instructions for use.

- Air in the infusion set can lead to an insufficient insulin delivery.

Follow the instructions in these instructions for use on how to fill the infusion set completely or how to remove any air bubbles.

Risk of hyperglycaemia or inflammation

Mechanical stress, for example, pressure or

force, on the infusion set or the infusion site can lead to an insufficient insulin delivery or inflammation.

Avoid mechanical stress on the infusion set and on the infusion site.

Risk of inflammation

In rare cases, the product can cause unexpected body reactions.

If you experience skin irritation or redness, contact your healthcare professional.

Risk of infection, inflammation, or scarring

If you use infusion set components for a longer period than intended, infection, inflammation, or scarring can occur.

Replace the cannula at least every 3 days.

Replace the transfer set at least every 6 days.

If your blood glucose level unexpectedly increases or if an occlusion alarm occurs, check the infusion

set for occlusion or leaks. If you are not sure, whether the infusion set is working properly, replace the infusion set immediately.

If your infusion site is inflamed, replace the infusion set immediately. Select a different infusion site during the healing process.

If the head set or the transfer set was dropped or exposed to other mechanical stress, check it for damages, for example, cracks. If it is damaged, do not use it.

Check the connection between the luer-lock connector (J) and the insulin pump adapter every 3 hours during the day and before you go to bed.

Check your blood glucose level at least 4 times a day.

Check at regular intervals that the adhesive pad does not peel off from the infusion site and that the cannula is fully inserted. If the cannula is not

fully inserted anymore, replace the cannula.

After you replace the infusion set or the cartridge and before you connect the transfer set to the cannula, check whether the transfer set is completely filled with insulin.

Suitable Infusion Sites

Select an infusion site with sufficient subcutaneous tissue, for example, on the abdomen or the back of your upper arm. See the illustration on page 2. The infusion site must be at least 5 cm away from the navel and the last infusion site used and must not be over a bone. The infusion site must not be exposed to rubbing or bumping, for example, through a waistband, belt, or safety belt. The infusion site must not have skin changes, for example, scars, injuries, bruises, rashes, or tattoos.

Rotate between different infusion sites each time

en

9

you replace the cannula to avoid tissue changes from insulin. Think of a pattern to select new infusion sites.

Attaching the Infusion Set

When you use the infusion set for the first time, your healthcare professional must be present.

- Wash your hands.
- Disinfect the chosen infusion site as instructed by your healthcare professional.
- Let your skin dry completely.

PRECAUTION

Risk of infection

If you do not disinfect the infusion site, an infection can occur.

Disinfect the selected infusion site as instructed by your healthcare professional before you insert the cannula into the skin.

10

See illustrations on page 1.

- 1** Hold the head set by the handling aid (A).
Pull off the needle cover (G).
Do not touch the introducer needle and the cannula to avoid injury and contamination.
- 2** Pull off both parts of the protective film (D) from the adhesive pad (C).
- 3** Form a skin fold at the chosen infusion site.
Insert the introducer needle quickly, in a single motion, and at a 90° angle into the skin fold.
- 4** Press the adhesive pad (C) firmly onto the skin.
Check that the adhesive pad is properly attached to the skin.
- 5** Hold the adhesive pad (C) in place.
Press the side wings (B) of the handling aid (A) together.

Pull off the handling aid (A) in the direction of the arrow.

- 6** Hold the cannula casing (H) in place.
Pull out the introducer needle straight upwards and discard it.

7 **PRECAUTION**

Risk of hypoglycaemia, hyperglycaemia, or infection

If the luer-lock connector is damaged, the insulin delivery can be insufficient or uncontrolled or germs can enter the infusion set.

To avoid damaging the luer-lock connector, do not use any tools to screw the luer-lock connector into the insulin pump adapter. Only screw by hand.

Screw the luer-lock connector (J) of the

transfer set tight into the insulin pump adapter.

• **WARNING**

Risk of hypoglycemia

If you fill the transfer set while the infusion set is connected to your body, an uncontrolled insulin amount can be delivered.

Never fill the transfer set while the infusion set is connected to your body.

Fill the transfer set as described in the User's Manual of the insulin pump until insulin drips from the needle inside the connector (L). You must fill up the transfer set completely before you connect it to the cannula. The transfer set

en

must not contain any air bubbles. Remove any air bubbles while filling the transfer set.

Fill volume of the transfer set (U100 insulin)

Length	Units	Volume
30 cm / 12"	approx. 6 U	approx. 60 µL
60 cm / 24"	approx. 11 U	approx. 110 µL
80 cm / 31"	approx. 14 U	approx. 140 µL
110 cm / 43"	approx. 18 U	approx. 180 µL

- 8** Hold the cannula casing (H) in place. Push the connector (L) onto the cannula casing (H). For a tight connection, make sure that you hear and feel the connector click into place.
- Fill the cannula casing (H) with a standard bolus of 1.0 U insulin (10 µL). Refer to the User's Manual of your insulin pump for how to

program a standard bolus.

If air is still in the infusion set, disconnect the transfer set from the cannula. Repeat the filling procedure for the transfer set as described in the User's Manual of the insulin pump. If the air is completely removed, reconnect the transfer set to the cannula.

Refer to the User's Manual of your insulin pump for how to start insulin delivery.

Disconnecting and Reconnecting the Insulin Pump

This infusion set allows you to disconnect your insulin pump, for example, for showering and swimming. You can reconnect the insulin pump without having to replace the infusion set.

Consult your healthcare professional on how to compensate for any missed insulin while you are

disconnected from the insulin pump.

Disconnecting the insulin pump

- Hold the adhesive pad (C) in place.
- Press the side wings (M) of the connector (L) of the transfer set together.
- Pull the connector (L) off the cannula casing (H).
- Push the disconnect cover (I) onto the cannula casing (H).

Reconnecting the insulin pump

- Hold the adhesive pad (C) in place.
- Press the side wings of the disconnect cover (I) together.
- Pull the disconnect cover (I) off the cannula casing (H).
- Push the connector (L) of the transfer set onto the cannula casing (H). You must hear and feel

the connector click into place.

When reconnecting the insulin pump, the cannula casing is still filled with insulin. Therefore, you do not need to fill the cannula casing with an insulin bolus.

Replacing the Cannula or the Transfer Set

- Put your pump into STOP mode.
- Disconnect the transfer set from the cannula casing (H), see section “Disconnecting the insulin pump”.

Replacing the cannula

- Carefully pull off the adhesive pad (C), working round the outside edge into the center.
- Pull the cannula (E) out of the skin and completely pull off the adhesive pad simultaneously.

en

- Insert a new cannula. Start with step **1** in the section “Attaching the Infusion Set”.

Replacing the transfer set

- Unscrew the luer-lock connector (J) of the transfer set from the insulin pump adapter.
- Do the steps **7** to **8** in the section “Attaching the Infusion Set” with the new transfer set.

The cannula casing is still filled with insulin. Therefore, you do not need to fill the cannula casing with an insulin bolus.

Discarding Used Parts of the Infusion Set

PRECAUTION

Risk of infection

A used product that has come into contact with human body fluids can transmit infections.

Discard all used parts of the infusion set as infectious material according to local regulations. For more information see the following instructions.

Discard the pointed introducer needle in such a way that it does not cause injury to anyone. For example, place the introducer needle into a sharps container.

For information on how to discard the product correctly, contact your local council or authority.

Operating and Storage Conditions

Use the infusion set between 5 and 40 °C.

Store the infusion set in the unopened blister at room temperature. Temporarily, the temperature can be between 5 and 45° C. Keep the product dry and away from sunlight.

Reporting of Serious Incidents

For a patient/user/third party in the European Union and in countries with identical regulatory regime; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

Customer Support

Australia

Roche Diabetes Care Australia Pty. Limited
Insulin Pump Support: 1800 633 457
australia.insulinpumps@roche.com
www.accu-chek.com.au

Roche Diabetes Care South Africa (Pty) Ltd.
Hertford Office Park, 90 Bekker Road
Vorna Valley, 1686

South Africa

Tel: +27 (11) 504 4600
Email: info@accu-chek.co.za
www.rochediabetescaremea.com

Singapore

Accu-Chek ExtraCare line: 6272 9200
www.accu-chek.com.sg

Distributed in the United Kingdom by:

Roche Diabetes Care Limited
Charles Avenue, Burgess Hill
West Sussex, RH15 9RY, **United Kingdom**

Accu-Chek Pump Careline ¹⁾:

UK Freephone number: 0800 731 22 91

ROI Freephone number: 1 800 88 23 51

¹⁾ calls may be recorded for training purposes

Some mobile operators may charge for calls
to these numbers.

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