ACCU-CHEK, ACCU-CHEK TENDERLINK and AKKY-CHEK
are trademarks of Roche.

© 2019 Roche Diabetes Care

Made in Mexico

Distributed by:

Roche Diabetes Care GmbH
Sandhofer Strasse 116
68305 Mannheim, Germany

Unomedical a/s
Aaholmvej 1–3, Østved
4320 Lejre, Denmark


ACCU-CHEK® TenderLink







Instructions for use

Infusion set

ACCU-CHEK®


Instructions for use	4	Kasutusjuhised	224
Gebrauchsanweisung	20	הוראות שימוש	255
Mode d'emploi	36	Інструкція з використання	256
Istruzioni per l'uso	52	Használati útmutató	272
Instrucciones de uso	68	Naudojimo instrukcijos	288
Gebruiksaanwijzing	84	Lietošanas instrukcija	304
Instruções de utilização	100	Instrucțiuni de utilizare	320
Návod k použití	116	Инструкции по использованию	336
Instrukcja obsługi	130	Návod na používanie	352
Brugsanvisning	146	Navodila za uporabo	368
Bruksanvisning	162	إرشادات الاستخدام	399
Bruksanvisning	178		
Käyttöohje	192		
Οδηγίες χρήσης	206		

Instructions for use	en
Gebrauchsanweisung	de
Mode d'emploi	fr
Istruzioni per l'uso	it
Instrucciones de uso	es
Gebruiksaanwijzing	nl
Instruções de utilização	pt
Návod k použití	cs
Instrukcja obsługi	pl
Brugsanvisning	da
Bruksanvisning	sv
Bruksanvisning	no
Käyttöohje	fi
Οδηγίες χρήσης	el


ACCU-CHEK® TenderLink

Infusion set

head set

transfer set

Intended use

The Accu-Chek TenderLink infusion set is intended for the subcutaneous infusion of insulin, administered with microdosage insulin pumps.

Package 1:

- | | |
|--|--|
| A protective cover | B soft cannula |
| C introducer needle | D side clips of the introducer needle |
| E self-adhesive | F front half of the adhesive backing |
| G back half of the adhesive backing | H cannula housing |
| I tubing | J connector with connector needle |
| K side clips of the connector | L circular protective cap |
| M disconnect cover | N luer-lock connector |

Package 2:

Components A–H + M from package 1.

Contents of package I

Package 1 only.

Contents of the cannula package for the infusion set

Package 2 only.

WARNING

Read the instructions carefully before use. It is recommended that you consult your insulin pump manual for information regarding pump therapy.

- When using the infusion set for the first time, your doctor or healthcare team must be present.
- This device remains sterile until the package is opened or damaged. Do not use if the sterile package has been opened or damaged.
- Protect the product from excessive humidity, sunlight, and heat. Store at room temperature.
- Discard your infusion set after use according to local regulations. Do not clean or re-sterilize.

- Reuse of the infusion set may cause infection, site irritation, or damage to the cannula/needle. A damaged cannula/needle may lead to inaccurate medication delivery.
- The infusion set is for subcutaneous use only.
- Always screw the luer-lock connector firmly into the pump adapter. If this is not done, a loose connection between the infusion set and the pump could cause leakage. Avoid using any tools to screw the luer-lock connector into place, otherwise damage could result to the luer-lock connector.
- Do not leave air bubbles in the infusion set. Prime thoroughly.
- A small percentage of the Accu-Chek TenderLink infusion sets using soft cannulas can accidentally become crimped during insertion or displaced while in use. If this should occur, use a new infusion set immediately.
- Check frequently to make sure the soft cannula remains firmly in place. Because the cannula is so soft, you will not notice any pain if it is pulled out. The soft cannula must always be completely inserted in order to receive the

- full amount of insulin. Replace the infusion set as soon as the self-adhesive loosens.
- Do not reinsert the introducer needle into the soft cannula. Reinsertion could cause tearing of the soft cannula which would result in unpredictable insulin flow.
- If you have to disconnect your infusion set, pay particular attention to hygiene.
- Never use a soft cannula for more than 72 hours (3 days). Never use the transfer set for more than 6 days. When changing the cartridge, make sure that the transfer set is completely primed before connecting to the head set.
- If the infusion site becomes inflamed, replace the infusion set and use a new site.
- Do not put disinfectants, perfumes, deodorants, or insect repellent on your infusion set as these may affect the integrity of the infusion set.

- Never prime the transfer set or attempt to free clogged tubing while the transfer set is connected to the head set. You may accidentally infuse an uncontrolled quantity of insulin.
- Avoid excessive stress on the infusion site and on the infusion set. Pay attention when carrying heavy weights.
- Before inserting the infusion set, prepare the infusion site in accordance with the instructions of your healthcare professional.
- Check your blood glucose level 1 to 3 hours after inserting your infusion set and inspect the infusion site on a regular basis. Do not insert a new infusion set just prior to bedtime. Check blood glucose frequently (check with your physician).
- If there is an unexplained increase in your blood glucose level, causing it to be too high, or if an occlusion alarm occurs, check for clogs and leaks. If in doubt, change your infusion set.

Insert the infusion set

Picture 1

Choose a site away from the waistline, bones, scar tissue, belly button and recent infusion sites. Wash your hands and clean the infusion site according to the instructions from your healthcare professional. Allow the infusion site to dry before the inserting the infusion set.

Remove the front half of the adhesive backing (F) from the self-adhesive (E) and pull the protective cover (A) off the introducer needle (C). Grasp the infusion set between the thumb and middle finger of one hand behind the side clips and pull back the edge of the self-adhesive with the index finger.

WARNING

Ensure clean application. Do not touch the soft cannula or the introducer needle.

Picture 2

With your free hand, pinch together a bit of subcutaneous tissue at the disinfected site. Insert the introducer needle (C) under the skin at an angle of 20–45°, proceeding quickly and in one motion. Watch the soft cannula (B) while inserting to assure that it goes in smoothly without kinking or crimping.

⚠ WARNING

People with very small or large amounts of subcutaneous fatty tissue should take care when choosing the angle of insertion or the cannula length. The cannula could be lodged in the muscular tissue or epidermis, where absorption of the insulin is limited or impossible. Ask your healthcare professional for more information.

Make sure that the soft cannula is fully inserted into the subcutaneous tissue and not kinked at the skin surface. Insert the introducer needle in a single, smooth movement, otherwise the soft cannula may kink or crimp upon insertion, causing reduced insulin delivery or a blockage.

Picture 3

Press the self-adhesive (E) firmly onto the skin in order to ensure secure seating of the infusion set. Place one finger in front of the inspection window on the infusion set to fix it. Press the side clips of the introducer needle (D) together gently and at the same time pull the introducer needle (C) out of the skin.

WARNING

Discard your introducer needle and infusion set after use according to local regulations.

Picture 4

Remove the back half of the adhesive backing (G) from the self-adhesive (E) and press the self-adhesive firmly onto the skin.

Picture 5

Attach the luer-lock connector (N) to the pump by screwing it into the pump adapter. Make sure there is a cartridge containing sufficient insulin in the cartridge compartment. Remove the circular protective cap (L) by pressing the side clips gently inwards. Follow the pump manufacturer's User's Manual to prime the transfer set until insulin drips from the connector needle tip. Do not leave air bubbles in either the cartridge or the transfer set.

Picture 6

Connect the transfer set to the newly inserted head set with the bulging side upwards. A click must be heard when the connection engages. Prime the empty space of the cannula housing (H) with a bolus of 0.7 units of U100 insulin (about 7µl). Start insulin delivery with the pump.

⚠ WARNING

Change the head set every 2 to 3 days (max. 72 hours), otherwise you risk insulin resistance or inflammation of the infusion site.

Never use the transfer set for more than 6 days.

You must prime the empty space of the newly inserted head set after connecting the transfer set. Failure to do so results in missed insulin.

Always check your blood glucose level 1 to 3 hours after introducing the set.

If the coupling does not click noticeably and audibly together, use a new head set and a new transfer set.

Do not use any tools to screw the luer-lock connector into place.

Priming volume:

30 cm/12"	approx 7 units U100 insulin	approx 70 µl
60 cm/24"	approx 10 units U100 insulin	approx 100 µl
80 cm/31"	approx 13 units U100 insulin	approx 130 µl
110 cm/43"	approx 18 units U100 insulin	approx 180 µl

Disconnect the infusion set

The infusion set allows you to temporarily disconnect from your pump for showering, swimming, etc., without changing your infusion set.

WARNING

Consult with your healthcare professional on how to compensate for any missed insulin while you are disconnected. Carefully monitor your blood glucose level while disconnected from the pump and after you reconnect.

Step 1

Stop the pump. Hold the self-adhesive (E) in place, gently press the side clips of the connector (K) together and pull the connector (J) from the cannula housing (H).

Step 2

Put the disconnect cover (M) on the cannula housing.

Reconnect the infusion set

WARNING

Before reconnecting, ensure that a drop of insulin appears at the tip of the connector needle.

Step 1

Prime the infusion set until insulin comes out of the connector needle.

Step 2

Hold the self-adhesive (E) in place and remove the disconnect cover (M) from the cannula housing (H).

Reattach the connector (J) to the cannula housing. An audible “click” ensures the connector is locked.

WARNING

It is not necessary to give an additional bolus.

United Kingdom

Distributed in the United Kingdom by:

Roche Diabetes Care Limited

Charles Avenue, Burgess Hill

West Sussex, RH15 9RY, **United Kingdom**

Accu-Chek Pump Careline ¹⁾:

UK Freephone number: 0800 731 22 91

ROI Freephone number: 1 800 88 23 51

¹⁾ calls may be recorded for training purposes

Some mobile operators may charge for calls to these numbers.

burgesshill.insulinpumps@roche.com

www.accu-chek.co.uk

www.accu-chek.ie

Australia

Roche Diabetes Care
Australia Pty. Limited
Pump Support: 1800 633 457
www.accu-chek.com.au

New Zealand

Accu-Chek Enquiry Line:
0800 80 22 99
Pump Hotline: 0800 696 696
www.accu-chek.co.nz


Roche Diabetes Care South Africa (Pty) Ltd.
Hertford Office Park, 90 Bekker Road
Vorna Valley, Midrand

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